

## Cloning of gRNAs into pLENTICRISPRv2-GFP plasmid

**Aim:** To generate pLENTICRISPRv2-GFP plasmids containing sgRNAs of interest for CRISPR Cas9 mediated knock out of a gene of interest.

### Method:

- For this protocol we will use the plasmid pLENTICRISPRv2-GFP
  - This plasmid was modified to express GFP by the Feldser lab using a plasmid originally from the Feng Zhang lab which is available to purchase from addgene (<https://www.addgene.org/82416/>)
  - This plasmid is useful as it will contain the sgRNA of interest, cas9 and a GFP marker all in one plasmid
  - When transfected into HEK293T cells along with packaging and envelope coding plasmids (psPax2 and pMD2.G), lentiviral particles will be produced to transduce hard to transfect cell lines
- The gRNA sequence can be designed with an online program or sgRNAs can be used from previously published sequences.
- For use in the pLENTICRISPRv2-GFP plasmid add CACCG to the 5' end of the forward primer and then create a complement of the gRNA sequence and add on a C to the beginning and CAAA to the end. Then do a reverse to obtain the oligo to be ordered
  - This is so the annealed oligos can be ligated into the linearized pLENTICRISPRv2-GFP plasmid (digested with Esp3i enzyme (also known as BsmBi enzyme), ThermoFisher Scientific # FD0454). As shown below

### Target Guide Sequence Cloning Protocol

In order to clone the target sequence into the lentiCRISPR backbone, synthesize two oligos of the form:

Target Sequence: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 NGG PAM

Oligo 1 → 5' - CACCGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN - 3'

3' - CNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNCAAA - 5' ← Oligo 2

*Example oligo design:* Note that the NGG PAM is **not** included in the designed oligos.

Target Sequence: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 NGG PAM

Genomic 5' - . . . GACCACAGTCTGATCAGTTTTTCCTTGGGCTGCAA . . . - 3'

Sequence 3' - . . . CTGGTGTCAGACTAGTCAAAAGGAAACCGACGTT . . . - 5'

Oligo 1 → 5' - CACCGCAGTCTGATCAGTTTTTCCTT - 3'

3' - CGTCAGACTAGTCAAAAGGAAACAAA - 5' ← Oligo 2

*Oligonucleotide ordering tips:* Standard de-salted oligos (usually the most inexpensive synthesis) are sufficient for cloning. If not already resuspended, dilute each oligo to 100µM in sterile water or TE.

\* \* \* \* \*

**Figure 1.** Taken from Zhang lab protocol for construction of pLENTICRISPR plasmids  
[https://media.addgene.org/data/plasmids/52/52961/52961-attachment\\_B3xTwla0bkYD.pdf](https://media.addgene.org/data/plasmids/52/52961/52961-attachment_B3xTwla0bkYD.pdf)

- For example the below oligos were ordered to knockout mouse PD-L1

Gene name	gRNA sequence	Oligonucleotide sequence	PMID
CD274 / PD-L1 Forward	AATCAACCAGAGAATTTCCG	CACCGAATCAACCAGAGAATTTCCG	28723893
CD274 / PD-L1 Reverse	TTAGTTGGTCTCTTAAAGGC	AAACCGGAAATTCTCTGGTTGATTC	

- I order oligos dry from ThermoFisher scientific and then reconstitute in nuclease free water at 100  $\mu$ M
- The below annealing and ligation protocols were adapted for this plasmid from researcher Dr Abdullah Khan, University of Birmingham (@abattacks). Thanks Ab!

### Annealing of oligos

- First you need to anneal the oligonucleotides, this can be done without phosphorylation
- In a PCR tube assemble the following:

Oligo 1 Forward (100 $\mu$ M)	2 $\mu$ L
Oligo 2 Reverse (100 $\mu$ M)	2 $\mu$ L
5 x Ligase buffer (T4 cat no# 15224017)	4 $\mu$ L
Nuclease free water	12 $\mu$ L
TOTAL	20 $\mu$ L

Then heated to 95°C for 5mins and then turn off the PCR machine or heat block to allow the samples to cool down gradually with the tubes inside

### Digestion and ligation one step reaction

If you trust your pipettes to manage these small quantities go ahead, if not just double up on all the quantities.

pLENTICRISPRv2-GFP	0.5 $\mu$ L (25 ng)
Annealed oligos (1:1 dilution)	0.5 $\mu$ L
Esp3i (BsmBi) (Thermo, cat no. FD0454)	0.25 $\mu$ L
Fast digest buffer 10x	0.5 $\mu$ L
T4 ligase (Thermo, cat no. 15224017)	0.25 $\mu$ L
DTT 20 mM	0.25 $\mu$ L
Nuclease free water	2.75 $\mu$ L
TOTAL	5 $\mu$ L

Incubate at 37 °C for 2 hours and then transform 2 µL into Stab13 competent *E.coli* (I used home made competent cells, but can be bought from ThermoFisher Scientific).

I also included negative controls with all the above components but no oligos.

I got plenty of colonies and 24 out of 25 were correct by colony PCR. So I just sequenced 1 clone each construct and the large majority were the correct clone.

**Notes:**

- For sequencing of this vector I used the below forward primer
  - o TACGTGACGTAGAAAGTA
- For colony PCR to ID clones with an inserted gRNA I used the above forward sequencing primer and also the reverse cloning primer of each construct.